

O-INT^{ELLIGENCE} SERIES SMART SERIES



OT SERIES OXYGEN CONCENTRATOR

Model	OT-Elite3	OT-Elite5	OT-Elite6	OT-Elite8	OT-Elite10
Parameters					
Operating System	Timed Cycle/Rapid Pressure Swing				
Power Consumption (W)	350			480	
Electrical Requirement (V/Hz)	220/110±22V 50/60±1				
Delivery Rate(L/min)	0-3	0-5	0-6	0-8	0-10
Oxygen Concentration (%)	93±3%				
Outlet Pressure (MPa)	0.05-0.08				
Sound Level (dBA)	≤40			≤50	
Operating Environment (°C)	10-40				
Humidity Range (%)	≤80				
Equipment Class and Type	Class II a				
Safety Alarms & Features	Power Failure, High Pressure / Low Pressure, Low Oxygen Concentration (Model OT-ELITE, OT-PRO), Temperature, Current overload or line surge shutdown, Thermal Switch, 40 psig Pressure Relief Valve				
PURITY ALARM (REAL TIME)	OT-Pro3/Pro5/Pro8/Pro10				





O-INT^{ELLIGENCE} SERIES

GREY WHITE

PRODUCT HIGHLIGHT

- Real Time Oxygen concentration display
- Multi-figures Display
- Large Color LCD screen
- Clinical Software for nursing home & doctors
- Build in 2G Micro SD card
- Low noise version
- Unique Maintenance Alarm design for 4000 hours working



O-INT^{ELLIGENCE} SERIES
SMART SERIES

